



VIT[®]

Vellore Institute of Technology

(Deemed to be University under section 3 of UGC Act, 1956)

School of Advanced Sciences

VIT, Vellore – 632 014, Tamil Nadu

Contact: 0416-220 2463 & ftir.sas@vit.ac.in

Requisition Form for FT-IR Analysis (Internal Users)

Name:

Registration No.:

Guide Name:

School:

Email Id:

Contact Number:

The number of samples:

Nature of sample: Solid/ Liquid/ Film/Gel

Mode of analysis: **KBr Disc / ATR Mode:**

Purpose: Project / Ph D work / SET / J component

Sample Codes:

Are they air and moisture sensitive? : **Yes/ No**

Signature of the Student

Signature of the Guide

User Instructions

- ✓ Before submitting the samples for FT-IR analysis, you should confirm your slot date and time
- ✓ Minimum requirement of samples 3-5 mg for fine powder, 0.5 ml for liquid, and 2cm*2cm for film